




MASS SPECTROMETRY  
LABORATORY



CHEMISTRY  
UNIVERSITY OF  
TORONTO



## AIMS INSTRUMENTATION & SAMPLE REPORT DOCUMENTATION

### JEOL AccuTOF GCx plus



**AIMS MASS SPECTROMETRY LABORATORY**  
DEPARTMENT OF CHEMISTRY, UNIVERSITY OF TORONTO  
80 ST. GEORGE STREET, TORONTO ON, M5S 3H6

## **JMS-T200 AccuTOF Gcx-plus**

Manufacturer: [JEOL USA Inc.](#)  
Peabody MA

### ***Features***

Ion Sources: Electron Ionization (EI)  
Field Desorption (FD)  
Field Ionization (FI)

Ion Polarity: Positive and negative

Mass Range: 7-10 000 Da

Resolution: >10 000 mass resolution @  $m/z$  609

Mass Accuracy: +/- ~1.5 mDa

GC System: Agilent 7890B



### **AccuTOF 4G DART Results Format**

#### ***EI-MS Spectrum***


The *msAxel* Spectrum Report consists of a profile mass spectrum showing  $m/z$  versus intensity. Spectral peaks are labeled with  $m/z$  values rounded to two decimal places. The standard data acquisition range is  $m/z$  50-850 Da but printed mass spectra are normally expanded to display an appropriate region of interest. EI mass spectra are normally acquired in positive ion mode and the analyte species of interest appear as radical molecular ions  $[M]^+$ .


#### ***High Resolution Mass Spectrometry (HRMS) Report***

The EI-HRMS report is generated using the *Elemental Composition* feature implemented in the JEOL *msAxel* software package. The report provides the measured  $m/z$  value and a list of possible molecular formulae and calculated exact  $m/z$  values. All of the values contained in the report correspond to the *ionic species of interest, despite the table heading*. All ionic formulae are reported as charged species accounting appropriately for the mass of the electron (0.00054 Da).


#### **Suggested reporting format for compound characterization:**

"HRMS (EI-TOF+):  $m/z$   $[M]^+$ . calc'd for  $C_{18}H_{15}NF_4^+$   $m/z$  321.11351, measured  $m/z$  321.1142"





Spectrum



Elemental Composition

Parameters

Elements Set 1:

Tolerance: ±5.00 mDa  
 Electron: Odd/Even  
 Charge: +1  
 DBE: -1.5 - 40.0

Symbol	C	H	N	O	F
Min	0	0	0	0	4
Max	50	100	10	10	4

Results

Mass	Intensity	Formula	Calculated Mass	Mass Difference [mDa]	Mass Difference [ppm]	DBE
321.11416	196160.25	C18 H15 N F4	321.11351	0.65	2.03	10.0
		C5 H13 N10 O2 F4	321.11536	-1.19	-3.72	2.5
		C7 H15 N7 O3 F4	321.11670	-2.54	-7.90	2.0
		C15 H17 O3 F4	321.11083	3.33	10.37	5.5
		C9 H17 N4 O4 F4	321.11804	-3.88	-12.08	1.5
		C13 H15 N3 O2 F4	321.10949	4.67	14.55	6.0