



MALDI-TOF/TOF mass spectrometry

An Introduction

Contact : ms.support.us@bruker.com
Phone: 978-663-3660 Ext 5

1



Bruker Daltonics



Bruker MALDI-TOF(/TOF) Portfolio:



- **microflex LT** : Bench-top systems, 60 Hz nitrogen laser,
- **autoflex speed**: Floor standing, various configuration available (lin, ref, TOF/TOF), 2KHz for MS, 200Hz for MSMS;
- **ultrafleXtreme**: Top-of-the-line instrument, 2000 Hz MS, 1000 Hz for MSMS,

Bruker Daltonics







MALDI

7

Bruker Daltonics

MALDI stands for ...



**Matrix Assisted Laser
Desorption / Ionization**

invented by

**Tanaka (Nobel prize awarded 2002),
Hillenkamp and Karas**

in the mid '80s

8

Bruker Daltonics



How does MALDI work?



- Required items:**
- a **sample** to be analyzed
 - a **matrix** substance
 - a **laser**

MALDI preparation / MS analysis:

- 1.) **Mix sample and matrix** solutions in suitable ratio.
- 2.) Put a sub μ l aliquot of this mixture on the **target plate**.
- 3.) Let the mixture **co-crystallize**.
- 4.) Insert the target plate in the MALDI mass spec (running under **high vacuum**).
- 5.) **Shoot with the laser** on the MALDI preparation to **generate ions** (**singly charged** predominantly).
- 6.) **Collect, analyze and detect** the resulting **ions**.



9

How does MALDI work?



Fundamental papers on the principle of MALDI:

M. Karas, D. Bachmann, F. Hillenkamp

Analytical Chemistry, 57, 2935-2939 (1985)

K. Tanaka, H. Waiki, Y. Ido, S. Akita, Y. Yoshida, T. Yoshida

Rapid Communications in Mass Spectrometry, 2, 151-153 (1988)

R. C. Beavis, B. Chait, K.G. Standing

Rapid Communications in Mass Spectrometry, 3, 233-237 (1989)

M. Karas, M. Glückmann, J. Schäfer

Journal of Mass Spectrometry, 35, 1-12, (2000)

R. Zenobi and R. Knochenmuss

Mass Spectrometry Reviews, 17, 337-366 (1998)

10


Bruker Daltonics

What role does the **MALDI matrix** play?



The matrix transfers the energy needed for ionization from the laser light to the sample molecules.

Positive ionization mode:



Formation of alternative adducts depends on the presence of respective cations (either being ubiquitously present or actively added – depending on type of sample):
[M+Na]⁺; [M+K]⁺; [M+Cu]⁺; [M+Li]⁺; [M+Ag]⁺

11


Bruker Daltonics

What role does the **MALDI matrix** play?



The matrix transfers the energy needed for ionization from the laser light to the sample molecules.


Negative ionization mode:



12

Bruker Daltonics


Types of laser commonly used in MALDI:



Nitrogen laser:

pro: well structured energy profile


contra: slow (maximum 60Hz)



Nd:YAG laser:

pro: fast (up to >1000Hz)

contra: Gaussian energy profile (non-structured)



Smartbeam/Smartbeam II (modified Nd:YAG laser):

pro: fast (up to 2000Hz)

pro: well structured energy profile

Reference: A. Holle, A. Haase, M. Kayser, J. Höndorf, *Journal of Mass Spectrometry*, 41, 705-716 (2006)

13

Bruker Daltonics

Commonly used MALDI matrix substances:

Peptides: 4-Hydroxy- α -cyanocinnamic acid (**HCCA**)

Proteins: 2,5-Dihydroxyacetophenone (**DHAP**)

Sinapinic acid (**SA**)

2,5-Dihydroxybenzoic acid (**DHB**)

Glycans: 2,5-Dihydroxybenzoic acid (**DHB**)

Nucleic acids: 3-Hydroxypicolinic acid (**HPA**)

2,4,6-Trihydroxyacetophenone (**THAP**)

Why different matrices for different types of sample?

It's all about

- the **amount of energy needed to ionize** a particular sample compound

(individual matrices show specific „energy threshold“)

- the **stability** of a particular sample compound

(too „hot“ matrix may lead to non-desired fragmentation of sample compounds)

14

Bruker Daltonics

BRUKER

MALDI-TOF (Time of Flight)

15

Bruker Daltonics

BRUKER

MALDI-TOF

Basic principle of MALDI-Time-of-Flight based mass analysis:

The diagram illustrates the basic principle of MALDI-TOF mass spectrometry. It is divided into two main regions, both maintained at a high vacuum of 10^{-7} mbar.

- Acceleration Region:** On the left, the MALDI Ion Source (IS1) generates ions. These ions are accelerated through a series of electrodes, as indicated by the lightning bolt and the label "Acceleration".
- Field free TOF analyzer region (drift tube):** The ions travel through a drift tube where there is no electric field. The time taken for ions to travel this distance is the "Time Of Flight (depending on m/z)", represented by a stopwatch icon.
- Detector:** The ions are finally detected by a detector on the right.
- Mass Spectrum:** The detector output is plotted as Intensity versus m/z . The plot shows two distinct peaks, one pink and one green, representing different ion species.

16

Bruker Daltonics

MALDI-TOF



$$E_{\text{pot}} = zeU$$

$$E_{\text{kin}} = 1/2mv^2$$

$$zeU = 1/2mv^2$$

$$v = \sqrt{2zeU/m}$$

$$t = L\sqrt{m/2zeU}$$

This defines the potential energy at which all ions start from the MALDI target.

This equation defines the kinetic energy of ions after acceleration into the flight tube.

In the process of ion acceleration, energy is preserved but turns from potential into kinetic energy.

Transforming this equation shows the dependency of velocity of moving ions on their m/z value.

With the linear velocity v being defined as $L_{\text{drift tube}}/t_{\text{flight}}$, the dependency of the ions' flight time t_{flight} from m/z becomes obvious.

To put the basic MALDI-TOF separation principle into simple words:
The larger its m/z , the slower an ion will fly, the longer the measured flight time will be.

17

Bruker Daltonics


MALDI-TOF: Resolution



Resolution R defines, how well two peaks are separated from each other:

$$R = m/\Delta m$$

Δm is usually derived from a mass peak's full width at half maximum (FWHM), as shown below for the exemplaric analysis of a mixture of two compounds X and Y:




18

Bruker Daltonics

MALDI-TOF: Resolution

Example spectrum containing 3 singly charged $[M+H]^+$ signals at m/z 922Da, 923Da and 924Da:




$$\begin{aligned} \text{Resolution} &= m/\Delta m \\ &= 922.00852/0.07264 \\ &= \mathbf{12693} \end{aligned}$$

The resolution here is obviously more than sufficient to achieve baseline separation of the three compounds contained in the mixture.

19

Bruker Daltonics

MALDI-TOF: Linear mode



20


Bruker Daltonics

MALDI-TOF: Linear mode



Observation:

Ions of identical mass arrive at the detector at slightly different time points, causing peak broadening (i.e. limiting the resolution)



21

Bruker Daltonics


MALDI-TOF: Linear mode



Limited in resolution due to **spatial spread** and **energy spread**


Spatial spread:

- initial movement of ions towards different directions
- ions are desorbed from different z-coordinates due to heterogeneity in size of matrix crystals



Initial energy (=speed) spread:

- heterogeneous secondary reactions (ion-ion; ion-neutral)



Reference:

W. Ens, Y. Mao, F. Mayer, K.G. Standing,
Rapid Communications In Mass Spectrometry, 5, 117-123 (1991)

22

Bruker Daltonics

MALDI-TOF How to improve resolution???



- ❑ Optimized **matrix preparation homogeneity**
(minimizing spatial spread)
- ❑ **Pulsed ion extraction** for efficient ion focusing
(minimizing initial energy spread)
- ❑ Further ion focusing by means of a **reflector** TOF
setup (minimizing remaining energy spread)






23

Bruker Daltonics

MALDI-TOF How to improve resolution???




- ❑ Optimized **matrix preparation homogeneity**



24

Bruker Daltonics



How can I understand how PIE works?

PIE stands for **Pulsed Ion Extraction**. This means that the ions are not extracted as classically in a continuous extraction field, but are created in a zero extraction field, and then extracted after a certain delay time in a pulsed fashion.



This technique has two advantages: 1. The process is softer, as the ions are not accelerated through the dense plume of matrix molecules right after desorption, and 2. by choosing the right parameters a time focusing can be achieved. This together results in better resolution and better sensitivity, when using the correct parameters. On the other hand, wrong parameters can decrease resolution and sensitivity. Therefore the acquisition software greater than 3.0 allows for saving complete parameter sets. It is recommended to operate the system using parameter sets. Normally during installation some relevant sets have been stored on your system.

Method File Name	IS1	IS2	Lens	Delay
RP_3147	25	21.70	9.0	70ns
RP_6000	25	21.65	9.0	130ns
RP_8000	25	21.70	9.0	170ns

Method File Name	IS1	IS2	Lens	Delay
RP_1000	19	16.90	9.0	60ns
RP_2000	19	16.85	9.0	80ns
RP_3000	19	16.90	9.0	210ns
RP_6000	19	16.60	9.0	230ns
RP_12K	19	16.35	9.0	360ns
RP_17K	19	16.15	9.0	420ns

• Example only,

Bruker Daltonics



MALDI-TOF

Linear vs. reflector mode



FAQ:

If MALDI-TOF performed in reflector mode gives so much better resolution - why then use linear mode at all???

Answer:

Linear mode is used whenever analytes are not stable enough to survive the energetic stress which is inherent to passing the reflector (ions are decelerated/re-accelerated in the reflector by a high kV electric field within nanoseconds!!!).


Especially larger sized molecules, e.g. intact proteins, show limited stability when passing the reflector field, and may undergo serious fragmentation, which results in either badly resolved spectra (peak fronting due to non-resolved fragments) and/or drastic loss in sensitivity (low mass fragments will miss the reflector detector).


29

Bruker Daltonics

ultrafleXtreme:

On-axis linear detector: Outstanding performance in HMW measurement





MALDI-TOF

Where do these isotope peaks originate from?

Most elements are found in nature in form of different so called isotopes. Isotopes of an element have nuclei with the same number of protons (the same atomic number) but different numbers of neutrons.

Hydrogen 1 proton	^1H	^2H	^3H
Helium 2 protons	^3He	^4He	
Lithium 3 protons	^6Li	^7Li	
	Proton:	Neutron:	

However, both, protons and neutrons contribute to the weight of an atom, which explains the difference in weight of isotopes.

32 Bruker Daltonics

MALDI-TOF



Where do these **isotope peaks** originate from?

Isotope	Mass	[%] Abundance
1-H	1.007825	99.985
2-H (Deuterium)	2.014000	0.015
12-C	12.00000	98.90
13-C	13.00336	1.10
14-N	14.00307	99.63
15-N	15.00011	0.37
16-O	15.99491	99.76
18-O	17.99916	0.20
19-F	18.99840	100
23-Na	22.98977	100
31-P	30.97376	100
32-S	31.97207	95.03
34-S	33.96787	4.22
35-Cl	34.96885	76.77
37-Cl	36.96590	31.98
39-K	38.96371	93.26
79-Br	78.91834	50.69
81-Br	80.91629	49.31

Elements that are found in nature in form of only one single isotope, are called **monoisotopic elements**.

33

Bruker Daltonics


MALDI-TOF




Where do these **isotope peaks** originate from?

For a given molecule, the individual isotopes of all the elements contained in it finally yield a characteristic intensity distribution of isotopic masses. This is shown below by means of the isotopically resolved mass spectra of 3 compounds being different in size:


Element composition: $C_{41}H_{69}N_{13}O_{14}S$
 Monoisotopic mass $[M+H]^+$: 1000.4880
 Average mass $[M+H]^+$: 1001.1409



Element composition: $C_{112}H_{164}N_{29}O_{34}S_2$
 Monoisotopic mass $[M+H]^+$: 2524.1510
 Average mass $[M+H]^+$: 2525.8196



Element composition: $C_{253}H_{363}N_{55}O_{75}S$
 Monoisotopic mass $[M+H]^+$: 5404.6075
 Average mass $[M+H]^+$: 5407.9984



The **monoisotopic mass** is the sum of the masses of all the atoms present in a molecule using the mass of the most abundant isotope for each element.

The **average mass** of a molecule is the sum of elemental masses using the average weighted over all stable isotopes of each element contained in the molecule.

34

Bruker Daltonics

Calibration error, ppm calculation;



The calibration error (mass deviation) can be calculated by using the theoretical (m_{theo}) and measured (m_{measured}) mass of the unknown substance :

$$Error_{ppm} = \frac{m_{\text{measured}} - m_{\text{theo}}}{m_{\text{theo}}} * 10^6$$

To calculate the average calibration error over four measurements (m_1 , m_2 , m_3 and m_4) following formula needs to be used

$$D_{\text{avg}} = \frac{|m_{\text{theo}} - m_1| + |m_{\text{theo}} - m_2| + |m_{\text{theo}} - m_3| + |m_{\text{theo}} - m_4|}{4 * m_{\text{theo}}} * 10^6$$


Bruker Daltonics

MALDI-TOF: Calibration



Precision: Variation of values obtained from repetitive measurements performed under identical conditions (*random error*)


Accuracy: Deviation of a measured value from the reference value (*systematic error*)



36

Bruker Daltonics

MALDI-TOF: Calibration



$$t_{of} = t_{delay} + t_{acc} + t_{drift} \quad (1)$$

t_{acc} can be calculated as followed. The force to an charge q in an electric field E is

$$F = E q = M a \quad (2)$$

37

Bruker Daltonics

MALDI-TOF: Calibration



Where a is the mass dependent acceleration. It is related to the time t_{acc} and the distance d between target and grounded acceleration plate.

$$d = a/2 t_{acc}^2 \quad (3)$$

Replacing in (2) the electric field $E = U/d$, the charge $q = z e$, the number of charges z times the elementary charge e , and from (3) $a = 2d/t_{acc}^2$ than the time-of-flight through the acceleration region is

$$t_{acc} = d \sqrt{(2m/Uze)} \quad (4)$$

38

Bruker Daltonics

MALDI-TOF: Calibration



For the calculation of t_{drift} we have to do the following considerations. All ions start from the target at the same potential energy

$$E_{\text{pot}} = qU = zeU \quad (5)$$

Where q is the charge of an ion which z number of charges e . e is elementary charge and U is the high voltage applied to the target. Because of the conservation of energy, after acceleration the potential energy is transferred into kinetic energy

$$E_{\text{kin}} = 1/2 m v^2 \quad (6)$$

From the constant movement in a field-free region we know that the speed v is the distance L divided by the time t_{drift} . Replacing v in (6) by L/t_{drift} we get (as already shown earlier)

$$t_{\text{drift}} = L \sqrt{(m/2zeU)} \quad (7)$$

39

Bruker Daltonics

MALDI-TOF: Calibration



Combining of equation (4) – (7) the total time-of-flight is

$$t_{\text{of}} = t_{\text{delay}} + t_{\text{acc}} + t_{\text{drift}} = t_{\text{delay}} + d \sqrt{(2m/Uze)} + L \sqrt{(m/2zeU)}$$

TOF equation!

$$t = C_0 + C_1 \sqrt{m}$$

(Linear Mode)

$$t = C_0 + C_1 \sqrt{m} + C_2 m$$

(Quadratic or cubic enhanced)

Advanced reading: Calibration: Quadratic, linear and one-point calibration;

The determination of C_0 and C_1 is carried out by means of two signals of known standard masses m_1 and m_2 .

40

Bruker Daltonics

MALDI-TOF: Calibration




By measuring the time-of-flight t_1 and t_2 of two known masses m_1 and m_2 in a mass spectrum, C_0 and C_1 can be determined:

$$t_1 = C_0 + C_1 \sqrt{m_1}$$

$$t_2 = C_0 + C_1 \sqrt{m_2}$$



C_0, C_1



For a sufficiently narrow mass range, a linear approximation of t_{of} may be sufficient in terms of calibration accuracy (2 known calibrant masses are required here). However, for the calibration of an extended mass range, a quadratic polynomial provides a much better description of the relationship between flight time and m/z . Minimum 4 known calibrant masses are required when applying a quadratic calibration function!


41

Bruker Daltonics


MALDI-TOF: Calibration strategies



Step 1) External calibration




- apply calibration fit to sample spectra obtained from near neighbour spots



- calibrants of known mass cover mass range of interest
- m/z vs. flight time is fitted using a polynomial of varying order (depending on size of mass range to be calibrated and number of available calibrant signals, resp.)

Step 2) Internal re-calibration (optionally)



- * denotes compounds of known identity/mass
842.509 Da (trypsin artefact)
2211.104 Da (trypsin artefact)

Internal re-calibration allows for

- optimum mass accuracy due to compensation of spot-to-spot heterogeneities that typically cause mass errors after external calibration

42

Bruker Daltonics

MALDI-TOF: Calibration strategies



When to use which calibration polynomial?

Compass for Flex series 1.3 or higher:

	Min. number of calibrant peaks	Calibration polynomial to be used
external calibration:	2	linear
	4	quadratic
	6	cubic enhanced
internal re-calibration: (external pre-calibration: quadratic)	1	linear correction
	4	quadratic
	6	cubic enhanced
internal re-calibration: (external pre-calibration: cubic enhanced)	1	linear correction
	4	cubic enhanced

Note: For optimum mass accuracy, calibrants in general have to cover the entire mass range that is to be calibrated. Extrapolation of calibration functions will always have a negative effect on the resulting mass accuracy.

43


Bruker Daltonics

Mass Accuracy – Digitization Rate

How many data points are required to see peak features?



Bovine Insulin collected with 3 different digitization rates



Comparison of Peptide Mass Accuracy from 0.5, 1 and 2 GHz data

	0.5 GHz	1 GHz	2 GHz
Bradykinin	29 ppm	4.3 ppm	2.8 ppm
Angiotensin I	14.8 ppm	2.0 ppm	1.7 ppm

44

Bruker Daltonics

MALDI-TOF: Typical applications

Gel-based proteomics:
Protein ID by peptide mass fingerprinting (PMF)

- reflectron MALDI-TOF
- mass range: 700 ... 4000Da

Biological state A

Biological state B

- 1) 2D gel separation of proteins
- 2) Comparative image analysis
- 3) Excision of regulated spots
- 4) Enzymatic digestion
- 5) MALDI-TOF-PMF
- 6) Protein ID by database search (MASCOT, Phenyx etc.)

45 Bruker Daltonics

MALDI-TOF: Typical applications

Gel-based proteomics:
Protein ID by peptide mass fingerprinting (PMF)

Protein: Serum albumin precursor - Homo sapiens (Human) ALBU_HUMAN

Intensity coverage: 74.8% (676381 ions) | Sequence coverage MS: 49.8% | Sequence coverage MS/MS: 60% | pI: 5.9 | MW (kDa): 71.3

140	170	180	190	200	210	220	230	240
RFDFEIQE	ENFPAALVILIA	FAQFLQQCFP	EDRYKLNKV	TEFARTCVAD	ESAEKIQSL	HTLFGQLKLT	VATLRETYQK	NADCCAPKQ
NEETFLK	KYTELLEPH	FTYFLLEFL	FAERYKAFT	ECQQAADKIA	CLLPRLELP	DEGKAASAKQ	BLKCAALQYF	GERAFKAR
RLKCAID	RANLAKICE	WQDISSEAK	ECCERFLAK	SNCIARYVND	EMPADLPALA	ADPVEKQVC	KNYLKAQVY	LOHLYLIL
YADYQK	FPFLVEKPK	LKQKELPK	QLREYVQAL	LLVRYTQVQ	QSTYPLVQK	SNLGGVQSK	CKRPRASAK	PCARELVY
DALEVAT	YVPEFNAET	FTFADLITL	SEKEREIQKQ	TALVELVEIK	FRATKEQLKI	VNDQTAAYK	KKCAKURET	CFAREGRAL

46 Bruker Daltonics


MALDI-TOF: Typical applications

Gel-based proteomics:
Protein ID by peptide mass fingerprinting (PMF)



Why is MALDI-TOF the method of choice here?

- + speed (1 gel spot = 3-15 seconds of measurement)
- + fully automatized measurement
- + REALLY high throughput
- + sensitivity (low amol range)
- + good mass accuracy and sufficiently high sequence coverage to yield highly significant results




47
Bruker Daltonics

MALDI-TOF: Typical applications

MALDI-TOFMS profiling of microorganisms:



- linear MALDI-TOF
- mass range: 2000 ... 20000Da



Mix with MALDI matrix (HCCA),
prepare onto a MALDI target plate

Select a colony

?




Generate MALDI-TOF
profile spectrum

Fingerprint search against
proprietary spectra library

Identified species





Bacillus globigii

48
Bruker Daltonics

MALDI-TOF: Typical applications



MALDI-TOFMS profiling of microorganisms:

- + analysis speed
(5...15s per sample)
- + minimum sample prep
- + fully automated
- + highly significant results

49

Bruker Daltonics



MALDI-TOF/TOF

50

Bruker Daltonics

MALDI-TOF/TOF:

Principal scheme:

Ion source 1
CID cell
PCIS
LIFT
Ion source 2
PLMS
Reflector
Reflector detector

Ion path in TOF1 region (linear TOF)
Ion path in TOF2 region (reflector TOF)

Ion source 1 = MALDI ion source
Ion source 2 = LIFT re-acceleration cell
PCIS = Timed ion gate
PLMS = Post LIFT metastable suppressor

51 Bruker Daltonics

MALDI-TOF/TOF:

Analysis of a mixture containing 3 compounds (green, red, blue) being different in mass:

Source 1
CID cell
TOF1 region
PCIS
 $v_1 < v_2 < v_3$

Molecular ions are separated in TOF1 according to their mass.
Part of the ions undergo fragmentation in TOF1 region induced by:

- metastable laser induced decay (LID)
- collision induced decay (CID)

Most important:
Fragments and their precursor ions travel at the same velocity.

52 Bruker Daltonics

MALDI-TOF/TOF tandem MS:

LID vs. CID



LID: Laser-Induced Dissociation

Most straightforward way to peptide backbone fragmentation (b,y-type ions).

Used for **protein identification** by means of peptide sequencing.

CID: Collision-Induced Dissociation (high energy)

Additional side chain cleavages. Higher relative intensity of internal fragments. Overall shift of average fragment size towards lower mass.

Used as an option in special applications, e.g.:

- *denovo sequencing* (enhanced immonium ions)
- *differentiation* of isobaric aminoacids *L and I* by respective side chain cleavages
- *detailed glycan analysis* (monomer linkage positions determined by products of specific crossing cleavages)

53


Bruker Daltonics

MALDI-TOF/TOF tandem MS:



Analysis of a mixture containing 3 compounds (**green, red, blue**) being different in mass:

Selecting the **red** precursor ion for fragment analysis by MS/MS:



PCIS = Timed Ion Gate
(calibrated for time/mass correlation)

54

Bruker Daltonics

MALDI-TOF/TOF tandem MS:
 PCIS = Timed Ion Gate

BRUKER
 Bradbourny-Niellson

Source 1

PCIS

small fringing fields

l_{eff}

Braker Daltonics

55

MALDI-TOF/TOF tandem MS:
 PCIS = Timed Ion Gate

BRUKER

Ions to be rejected
 Ions to be selected
 Ions to be rejected

Braker Daltonics

56

MALDI-TOF/TOF tandem MS:
 PCIS = Timed Ion Gate

Ions to be rejected
 Ions to be selected
 Ions to be rejected

57 Bruker Daltonics


MALDI-TOF/TOF tandem MS:
 PCIS = Timed Ion Selector

Ions to be rejected
 Ions to be selected
 Ions to be rejected

58 Bruker Daltonics

MALDI-TOF/TOF tandem MS:

PCIS = Timed Ion Selector



Ions to be rejected
Ions to be selected
Ions to be rejected


59

Bruker Daltonics

MALDI-TOF/TOF tandem MS:

Analysis of a mixture containing 3 compounds (**green, red, blue**) being different in mass:

Red precursor ion selected for fragment analysis by MS/MS:



Source 1

Ion gate



LIFT
Ion Source 2



Rejected



Ion Source 2 (LIFT):
Re-acceleration of fragments and remaining molecular ions.
Re-focusing (resolution!!!) of the ions.

60

Bruker Daltonics








MALDI-TOF/TOF tandem MS:



Analysis of a mixture containing 3 compounds (green, red, blue) being different in mass:



PLMS:

Post-LIFT Metastable Suppressor deflects remaining non-fragmented precursor ions to avoid undesired metastable decay in TOF2 region, which would yield ghost fragment peaks detected at wrong mass.


67

Bruker Daltonics

MALDI-TOF/TOF tandem MS:





Analysis of a mixture containing 3 compounds (green, red, blue) being different in mass:



68

Bruker Daltonics



MALDI-TOF/TOF: Typical applications

Gel-based proteomics:

Protein ID based on PMF + data dependent MS/MS

Advantages of MALDI-TOF/TOF in gel-based proteomics:


- + speed (PMF + PFF analysis typically < 2min measuring time per gel spot)
- + amol sensitivity
- + automated, data dependent MS/MS workflow
- + improved significance of combined (PMF+PFF) datasets compared to PMF only data

71

Bruker Daltonics

MALDI-TOF/TOF: Typical applications

Analysis of posttranslational modifications: Phosphorylation




72

Bruker Daltonics

MALDI-TOF/TOF: Typical applications

Analysis of posttranslational modifications: Phosphorylation




73

Bruker Daltonics

MALDI-TOF/TOF: Typical applications

Analysis of posttranslational modifications: Phosphorylation




74

Bruker Daltonics

MALDI-TOF/TOF: Typical applications

LC-based proteomics: LC-MALDI analysis of complex samples

- 1) Optional pre-fractionation on protein level (gel; LC; IP etc.)
- 2) Enzymatic digestion
- 3) nanoLC of resulting peptide mixture
- 4) LC fraction collection on MALDI plate
- 5) offline MALDI-TOF/TOF analysis



75 Bruker Daltonics

MALDI-TOF/TOF: Typical applications

LC-based proteomics: LC-MALDI analysis of complex samples

- + reduced ion suppression due to separation upfront to the mass spectrometer
- + no time constraints during data acquisition (no time pressure caused by LC speed)
- + increased depth of analysis (more protein IDs from complex mixtures; increased sequence coverage for individual proteins)
- + data dependent, non-redundant MS/MS precursor selection (esp. useful for non-isobaric SILE experiments)
- + LC separation is „frozen“ on the MALDI plate and, thus, allows archiving of LC separated samples for later re-analysis
- a single LC-MS/MS experiment, when compared to online-ESI, is rather time consuming

76 Bruker Daltonics

BRUKER

MALDI-TDS:

Top-down sequencing of intact proteins
by means of MALDI-TOF/TOF

77


Bruker Daltonics

BRUKER

MALDI-TDS:


Top-down sequencing of intact proteins

N T E R M S E Q U E N C E C T E R M
 + MALDI matrix (1,5-DAN, sDHB)



BRUKER MALDI-TOF



In-Source Decay (ISD) upon MALDI process:







1804.156	1735.861	1512.527	1312.179	1162.865	1021.911	891.974	784.40	684.188	594.186	513.212	439.216	361.242	294.431	237.235	187.129	104.160
----------	----------	----------	----------	----------	----------	---------	--------	---------	---------	---------	---------	---------	---------	---------	---------	---------



78


Bruker Daltonics











87

Bruker Daltonics

Bruker Corporation Today

Life Science Analysis (LSA)



Bruker AXS —● Elemental Analysis

- X-ray Diffraction
- X-ray Crystallography
- X-ray Fluorescence
- EDS Microanalysis
- Optical Spectroscopy
- Combustion Analysis

Bruker BioSpin —● Magnetic Resonance

- NMR
- MRI
- EPR

Bruker Daltonics —● Mass Spectrometry

- LC MALDI-TOF(/TOF)
- LC Ion Trap MSⁿ
- GC and GC-(Qq) MS
- LC ESI-(Qq)-TOF, FTMS
- IMS
- ICP-MS

Bruker Optics —● Vibrational Spectroscopy

- FT-IR
- FT-NIR
- Raman
- TD-NMR

PRIME

The Best of LC CSI-MS and LC MALDI-MS





CaptiveSpray



amazon SL



amazon Speed ETD



UltrafleXtreme



maXis Impact



micrOTOF



amazon SL



Autoflex Speed TOF/TOF



micrOTOF QII



Autoflex Speed LRF



Autoflex Speed L



maXis 4G



Microflex LT



Microflex LRF




Proteineer FC



proteinscape



Threshold irradiance



Useful Analytical range

Resolution, mass centroid, and intensity of a peak as a function of laser power.

90 **Bruker Daltonics**

BRUKER

Reflector detector

The new reflector detector (auto-detector II) uses an electron multiplier instead of a multi channel plate (MCP). This electron multiplier has a shorter pulse width which allows to achieve higher resolution.

The detector can be tilted on two axes (tip and turn). To control the stepper motors the board UTLMC (attached to the lift plate) was replaced by the board UTSCB.

Because the electron multipliers need more operating voltage (up to 4000 V) than the MCPs (and also no gating is possible) the DSG module is replaced by two -5 kV power supplies.

Linear detector

The new linear detector uses also an electron multiplier but with less resolution than the previous MCP. The advantage here is the higher sensitivity and no saturation effects.








Learn more about differences between electron multiplier and MCP

Bruker Daltonics


BRUKER

Traditional laser beam path



Laser Attenuator Focus Ion Source Target

smartbeam™ laser beam path



Laser Attenuator Modulator Focus lens Pinhole Ion Source

"Slide"

92

Bruker Daltonics